

SLD TOILET LIFT SEAT

CLASSIC MODEL

Interactive
demonstration

Bedpan
accessory
(Optional)



Curve Lifting



Anti-injury seat



One button
control



SLD-TL-18-A

CLASSIC MODEL

Color: white

Material: ABS

Size: L60.6 * W52.5 * H71(cm)

Raising height: Front: 59(cm)

Rear : 80(cm)

Raising angle: 0-33°(max)

Waterproof level: IP44

Material: ABS

Net weight: 18kg

Function: lift up and down

Seat capacity: 200kg

Armrest loading capacity: 100kg

Working voltage: 110V~240V

Packaged size: L68 * W60 * H57(CM)



SLD TECHNOLOGY

SLD TOILET LIFT SEAT

BASIC MODEL

Interactive
demonstration

Lithium battery



Bedpan
accessory
(Optional)



One button
control



Anti-injury seat



Curve Lifting



SLD-TL-18-A

BASIC MODEL

Color: white

Material: ABS

Size: L61.6 * W52.5 * H79(cm)

Raising height: Front: 59(cm)

Rear : 80(cm)

Raising angle: 0-33°(max)

Waterproof level: IP44

Net weight: 18kg

Function: lift up and down

Seat capacity: 200kg

Armrest loading capacity: 100kg

Input Voltage voltage: 110v~240v

Working voltage: 24v DC

Packaged size: L68 * W60 * H57(CM)

SLD TOILET LIFT SEAT

COMFORTABLE MODEL

Backrest



Lithium battery



Anti-injury seat



Bedpan
accessory
(Optional)



Curve Lifting



One button
control



Interactive
demonstration



SLD-TL-18-A

COMFORTABLE MODEL

Color: white Material: ABS

Size: L61.6 * W52.5 * H79(cm)

Raising height: Front: 59(cm)

Rear : 80(cm)

Raising angle: 0-33°(max)

Waterproof level: IP44

Packaged size: L68 * W60 * H57(CM)

Net weight: 18kg

Function: lift up and down/backrest/
armrest

Seat capacity: 200kg

Armrest loading capacity: 100kg

Input Voltage: 110v~240v

Working voltage: 24v lithium battery



SLD TOILET LIFT SEAT

REMOTE CONTROL MODEL

Interactive
demonstration

Remote control



Lithium battery



Backrest



Bedpan
accessory
(Optional)



Anti-injury seat



One button
control



Curve Lifting



SLD-TL-18-A

REMOTE CONTROL MODEL

Color: white Material: ABS
Size: L61.6 * W52.5 * H79(cm)
Raising height: Front: 59(cm)
Rear : 80(cm)
Raising angle: 0-33°(max)
Waterproof level: IP44
Packaged size: L68 * W60 * H57(CM)

Net weight: 20kg
Function: lift up and down/backrest/
remote control(50m)
Seat capacity: 200kg
Armrest loading capacity: 100kg
Input Voltage: 110v~240v
Working voltage: 24v lithium battery

Suzhou SLD Mechanical Technology Co, .Ltd

No.	components name	project	Product type							Special comments
			Classic	Basic	Comfortable	Advanced	Remote control	Luxury	Smart washing and drying	
1	Toilet seat	Configuration	○	○	○	○	○	○	○	
2	Handrail	Configuration	○	○	○	○	○	○	○	
3	Battery	Configuration	—	○	○	○	○	○	—	
4	Arm rest	Configuration	—	—	○	○	○	○	○	
5	Backrest	Configuration	—	—	○	○	○	○	—	
6	One button call for urgent help	Configuration	—	—	—	○	—	○	Optional	
7	Remote control device	Configuration	—	—	—	—	○	○	—	
8	Washing and drying function	Configuration	—	—	—	—	—	—	○	
9	Roller for medical	Configuration	Optional	Optional	Optional	Optional	Optional	Optional	—	
10	Bedpan bracket	Configuration	Optional	Optional	Optional	Optional	Optional	Optional	—	
11	Crouch bracket	Configuration	Optional	Optional	Optional	Optional	Optional	Optional	—	

Special note: "○"means standard equipped, "—"means no equipped with this function,"Optional" Means the function is alternative equipped according to customer actually need



SLD TOILET LIFT SEAT

LUXURY MODEL

Remote control



Curve Lifting



Lithium battery



Backrest



One button control
(Optional)



Anti-injury seat



One button control



Armrest SOS alarm



Interactive demonstration



SLD-TL-18-A

LUXURY MODEL

Color: white

Material: ABS

Size: L61.6 * W52.5 * H79(cm)

Raising height: Front: 59(cm)

Rear : 80(cm)

Raising angle: 0-33°(max)

Waterproof level: IP44

Packaged size: L68 * W60 * H57(CM)

Net weight: 20kg

Function: lift up and down/backrest/
armrest SOS alarm/remote
control(50m)

Seat capacity: 200kg

Armrest loading capacity: 100kg

Input Voltage: 110v~240v

Working voltage: 24v lithium battery



SLD TOILET LIFT SEAT

INTELLIGENT WASH&DRY MODEL

Interactive
demonstration



SLD-TL-18-A

INTELLIGENT WASH&DRY MODEL

Seat heating



Curve Lifting



Massage
rinse



One button
control



Knob switch



Armrest SOS
alarm



Warm Air



Color: white

Material: ABS

Size: L61.6 * W52.5 * H79(cm)

Raising height: Front: 59(cm)

Rear : 80(cm)

Raising angle: 0-33°(max)

Waterproof level: IP44

Packaged size: L68 * W65 * H57(CM)

Net weight: 20kg

Function: lift up and down/wash/dry
clean hip / warming seat
night light

Seat capacity: 100kg

Armrest loading capacity: 100kg


Input Voltage: 110v~240v

Working voltage: 24v lithium battery



SLD TOILET LIFT SEAT

From the beginning SLD keeping creat comfortable elder's life by high technology



SLD TOILET LIFT SEAT

SENIOR MODEL

Interactive demonstration

Backrest



Lithium battery



Curve Lifting



Bedpan
accessory
(Optional)



Anti-injury seat



One button
control



Armrest SOS
alarm



SLD-TL-18-A

SENIOR MODEL

Color: white Material: ABS
Size: L61.6 * W52.5 * H79(cm)
Raising height: Front: 59(cm)
Rear : 80(cm)
Raising angle: 0-33°(max)
Waterproof level: IP44
Packaged size: L68 * W60 * H57(CM)

Net weight: 20kg
Function: lift up and down/backrest/
armrest SOS alarm
Seat capacity: 200kg
Armrest loading capacity: 100kg
Input Voltage: 110v~240v
Working voltage: 24v lithium battery