

EXCLUSIVE
NEWS

FOLDING +
ROTATION




FOLDING AND ROTATING TRANSFER BENCH

It is made entirely of aluminum and painted in epoxy resin with an elegant finish that adapts to any interior. It has a very reduced thickness of 22 mm which allows the installation on all cars and allows the door to be closed with the installed device reclosed and rotated. It is available for the driver and passenger side.

The new version with "rotation and retraction" is available, which frees the passage of the hand to the seat adjustment controls.

The fold-away folding transfer bench complies with the law nr. 104 of 5 February 1992 (Framework law for assistance, social integration and the rights of disabled people), therefore the right to tax relief is provided.
The RE-SEAT-02 folding transfer bench is adapted to the new regulation.



Transfer bench

Watch the video

